## Palajet®

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GB Working Instructions

Heraeus Kulzer GmbH  
Grüner Weg 11  
63450 Hanau (Germany)

info.lab@heraeus.com  
www.heraeus-dental.com

IM 24

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07.2011 Erstausgabe

### CONTENTS

<b>1. Designated Use</b> .....	13
<b>2. Equipment Characteristics</b> .....	13
<b>3. Declaration of Conformity</b> .....	14
<b>4. Starting the apparatus</b> .....	15
4.1 Installation .....	15
4.2 Preparation .....	15
4.3 Investment.....	15
4.4 Isolation .....	16
4.5 Preparing acrylic .....	16
4.6 Injecting .....	17
<b>5. Polymerising</b> .....	17
<b>6. Deflasking</b> .....	18
<b>7. Preparation for following processes</b> .....	18
<b>8. Cleaning – Maintenance – Repairs – Disposal</b> .....	18
<b>9. Safety Notes</b> .....	19
<b>10. Technical data</b> .....	20
<b>11. Contents of delivery and accessories</b> .....	20
<b>12. Name plate</b> .....	21
<b>13. Contacts in the countries</b> .....	22
<b>14. Document history</b> .....	22

## Validity of this document

This operation instruction manual is valid for Heraeus Kulzer Palajet – injection unit; Ref. 66020450; since serial-no. 600-07-2011; July 2011

## 1. Designated Use

The Palajet is a pneumatic injection apparatus for the production of dentures with PalaXpress® ultra and PalaXpress® cold-curing polymer or heat-curing polymers such as Paladon® 65.

Area of application: dental technology, such as laboratories, schools, universities.

The manufacturer does not accept liability resulting from use that is not in accordance with the instructions for use or unauthorised changes.

The apparatus is designed for continuous operation. Personal protective equipment, such as safety gloves, safety glasses etc., must be used at all times for protection from contamination with the processed materials.

## 2. Equipment Characteristics

- The Palajet is a pneumatic injection apparatus for injecting denture material into denture flasks.
- The apparatus is factory-set to a working pressure of 4 bar, which can be controlled with the pressure meter (see „Sources of error“).
- The pneumatic feed pressure of max. 10 bar must not be exceeded.
- The plunger speed during injection is 8-10 mm/s.
- The Palajet is designed and manufactured according to the most common safety regulations and guidelines.

### Operating Elements

1. Flask holder knob
2. Pressure module
3. Pressure plunger and cylinder
4. Pressure meter
5. Vent seal
6. Switch for operating the pressure plunger

## Description of Operation

During the injection process the mixed acrylic dough is pressed through the injection cylinder into the denture flask cavity by the pressure plunger. After the filling process is completed, the remaining denture acrylic will escape through the ventilation channel.

The ventilation channel is then sealed by the vent seal. (Please watch for deviations when using heat-curing polymers). The injection process runs automatically after setting the switch. We recommend the Palamat® elite polymerisation apparatus for polymerising the dentures.

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## 4. Starting the apparatus

On receipt of the apparatus, please check it for transport damage and report any damage to the transport company no later than 24 hours after delivery. Do not operate damaged apparatus under any circumstances.

### 4.1 Installation

The apparatus must be installed in dry areas only and on a solid foundation so that it cannot be tipped. Weight capacity of foundation: min. 20 kg; area: approx. 300 x 280 mm.

Before connecting the compressed air, check that the primary pressure is not more than 10 bar and not less than 4 bar. Also ensure that the compressed air does not contain oil or water. The apparatus is connected to the air supply with a hose (internal  $\varnothing$  8 mm).

### 4.2 Preparation

After connecting the compressed air, flip the switch down so that the pressure plunger releases the fill cylinder (pressure plunger moves down). Move the pressure module to the rear as far as it will go. Remove the fill cylinder from the apparatus by screwing anti-clockwise and remove it. Unscrew the fill cylinder cover from the fill cylinder.

### 4.3 Investment

Before embedding, make sure that both flask halves fit exactly (check for the same identification number). Lubricate the inside with Vaseline so that the plaster can later be easily removed. The Palajet singleflask (1 model) or the Palajet duoflask (1–2 models) can be used depending on the number of models. The spacer (placeholder) can be used to save plaster if the Palajet duoflask is used with only one model. If two models are used in the duoflask, the dentures must be connected with an injection channel on the opposite side of the injection opening in the tuber area.

The fabricated wax denture is embedded with the injection channel positioned at the maxillary tuberculum on the upper jaw or the retromolar triangle on the lower jaw. Remove excess plaster and smooth the surface. Clean excess plaster from the edges of the flask. Use class III plaster for embedding. After the plaster has hardened, attach the wax for the injection channel ( $\varnothing$  7 mm PalaXpress channel wax) and the ventilation channel (not required for Paladon 65).

#### **Applying channel wax:**

Position the ventilation channel (with  $\varnothing$  3 mm PalaXpress channel wax) in the area of the maxillary tuberculum or retromolar triangle (when using the single flask opposite the injection channel), choosing the shortest distance to the ventilation channel. Use a  $\varnothing$  7 mm connection channel with double flask or duoflask)

Insulate the plaster surface with Aislar®. Place the flask halves on top of each other and firmly screw them together with the included spanner (do not use a hammer on brass flasks). Fill in plaster through the flask opening so the teeth are covered and leave to harden. Remove excess material so the flask fits between the centring disks of the Palajet®.

#### 4.4 Isolation

After the plaster has hardened, preheat the flask in hot water (approximately 90°C (194°F)) without chemical additives in order to soften the wax, separate the flask halves after about 5 minutes and remove the remaining wax with clean hot water. Thinly insulate the hot plaster at the contact area to the denture plastic twice (allowing each layer to dry) with Aislar® (prevent formation of pools).

After isolation of the plaster roughen the basal area of the teeth with a coarse diamond (remove sanding dust). To improve the bond between the teeth and the denture plastic, apply Palabond® adhesive twice with a brush and allow each layer to dry for 30 seconds. After the second layer the adhesive will remain active for 10 minutes.

Close both flask halves firmly and insert the flask into the Palajet®. The sealing surface of the vent seal of the Palajet apparatus should not contact the opening of the flask's ventilation closure.

#### 4.5 Preparing acrylic

Prior to the mixing procedure unscrew the cylinder cover from the fill cylinder and insert the fill cylinder cartridge. In order to receive the acrylic dough, the seal disk with the sealing lip showing upwards must be pressed into the fill cylinder cartridge (as far as it will go). The cover cartridge will be inserted into the cylinder cover.

**Table 1**

Product	Mixing ratio
PalaXpress ultra	2 g : 1 ml 30 g : 15 ml
PalaXpress	2 g : 1 ml 30 g : 15 ml
Paladon 65	5 g : 2 ml 30 g : 12 ml

Pour the powder quickly into the reserved liquid, mix well for 30 seconds. Transfer acrylic dough to the filling cylinder immediately and cover.

Remove air bubbles by swinging the mixing cup. When working with PalaXpress ultra or PalaXpress pour the acrylic dough in a thin jet without bubbles into the prepared fill cylinder cartridge.

**Note**

To avoid inhomogeneities do not wipe off the mixing cup on the surface of the fill cylinder cartridge.

## 4.6 Injecting

### **PalaXpress/PalaXpress ultra:**

After filling the fill cylinder cartridge wait for the moment of injection. This moment is reached when the dough maintains a dull surface.

### **Paladon 65**

From the start of the mixing process the waiting period before injection at 23°C (73°F) is 20 minutes.

Waiting times are dependent on the room temperature and the mixed volume.

Close the fill cylinder with the cylinder cover and turn it into the Palajet. Lock the flask in position with the flask holder knob.

After pulling the pressure module as far as it will go, the switch is reversed while the fill cylinder is inserted and the denture acrylic is injected with help of the pressure plunger.

### **PalaXpress/PalaXpress ultra:**

As soon as material is pressed out of the ventilation channel, close the vent seal. This ensures that the pressure is built up in the flask. Leave the flask in the Palajet under pressure for 5 minutes (switch up). Set the timer to 5 minutes.

### **Paladon 65**

Leave the flask in the Palajet under pressure for 25 minutes (switch up). Set the timer to 25 minutes.

For operation of the timer see the separate instructions.

## 5. Polymerising

For the polymerisation of the material the flask must be removed from the injection apparatus. Place the switch in the starting position (down); the pressure plunger retracts from the injection cylinder. Once it is completely down, push the pressure module to the rear as far as it will go.

After unscrewing the fill cylinder the flask can be removed, placed on the flask filling basket and then inserted into the Palamat elite polymerisation unit. Disassemble the filling cylinder immediately and clean it. Acrylic residues can be removed much more easily from the filling cylinder and cover inserts after they have completely cured.

See „Preparation for following processes“.

### **PalaXpress/PalaXpress ultra:**

The time required for polymerisation is 30 min at 55°C (131°F).

**Paladon 65:**

Rapid polymerisation takes 30 min at 90°C (194°F).

Polymerisation in stages takes 2 h at 70°C–90°C (158°F–194°F).

Slow polymerisation takes 10 h at 90°C (194°F).

## 6. Deflasking

After polymerisation, the flask is removed together with the flask filling basket from the Palamat elite. Prior to deflasking the denture allow the flask to cool to room temperature, then deflask the denture. Remove the plaster from the flask halves with a plastic hammer.

**Note**

Do not use a metal hammer! Metal tools will deform the flask; this can result in inaccurate fitting. The final occlusion must be checked prior to removing the denture from the model!

## 7. Preparation for following processes

Prepare flask before reusing it and if necessary lubricate with Vaseline. Clean the sealing surfaces as well as the injection and ventilation channel from plaster and plastic remains in order to avoid raising the bite. Turn back the vent seal and remove any acrylic remains.

The fill cylinder and cover cartridge are constructed for multiple use (approximately 5 times). For better sealing and a safe pressure application a new sealing disk must be used for each injection procedure.

As noted above, acrylic residues can be removed from the filling cylinder and cover inserts much more easily after they have completely cured.

Take care that the cartridges are not damaged!

## 8. Cleaning – Maintenance – Repairs – Disposal

Clean the surface of the apparatus with a moist cloth. Make sure that no water can enter. Do not use aggressive cleaning agents. Check the water separator weekly and remove accumulated water as required. Check the compressed-air supply and the compressor if leaks occur. The apparatus may be damaged by moisture.

Inspection of the pneumatic components is recommended every two years (see „Service points“).

Use only approved original spare parts and accessories. The use of parts from other sources involves unknown risks and is not permitted under any circumstances.

The function and safety of the apparatus can only be guaranteed if all required tests, services and repairs are conducted by Heraeus service contractors or by technicians trained by the manufacturer.

Heraeus Kulzer GmbH will not accept liability for damages resulting from damage to the apparatus resulting from incorrect repair work that was not conducted by Heraeus service contractors or by technicians trained by the manufacturer, or if original spare parts or accessories were not used to replace parts or accessories.

The apparatus is designed for a service life of 10 years. In Germany please contact Heraeus Kulzer GmbH in Hanau for disposal of spare parts, components or the apparatus. In all other countries please contact the branch in your country or the authorised dealer.

**Sources of error**

Error sources	Remedy
Incorrect operating Pressure	<ul style="list-style-type: none"> <li>– Pull setting knob on the regulator- / filter unit upwards until raster snaps out</li> <li>– Adjust regulation of operating pressure to 4 bar and press setting knob downwards until it snaps in</li> </ul>
Too fast or too slow injection speed	<ul style="list-style-type: none"> <li>– Remove protective stirrup above the throttling check valve</li> <li>– Remove locknut on the throttling check valve</li> <li>– Adjust knurled head screw according to the desired speed and fasten the locknut</li> <li>– Replace protective stirrup</li> <li>– Hoisting speed = upper throttling check valve</li> <li>– Sink speed = lower throttling check valve</li> </ul>

**9. Safety Notes**

The safety of persons, operation and processing quality with this laboratory apparatus depends primarily on the actions of the persons who work with the apparatus.

Read the operating instructions carefully before using the apparatus and observe the directions for avoiding errors and resulting damage, particularly damage to health.

In addition to the directions in the operating instructions, the applicable national laws, regulations and guidelines must also be observed for installation and operation.

Maximum lift speed 50 mm/s. Flip the switch upwards only when the fill cylinder is inserted. Caution: danger of injury! Before working on the pressure module (pressure plunger) disconnect the apparatus from the compressed-air supply!

## 10. Technical data

Table 2

Pressure connection	max. 10 bar
Nominal pressure	4 bar
Plunger speed	0–50 mm/s
Compression force on the fill cylinder	ca. 200 kp
Height	ca. 400 mm
Width	ca. 290 mm
Depth	ca. 265 mm
Weight	13,5 kg

We reserve the right to make technical modifications.

**Note: Keep the operating instructions in an accessible position near the apparatus!**

Palajet® is a registered trademark of Heraeus Kulzer GmbH, 63450 Hanau

## 11. Contents of delivery and accessories

### Palajet base unit

- Injection cylinder including cover
- 1 flask (Palajet duoflask)
- 12 injection cylinder cartridges
- 12 cover cartridges
- 50 sealing disks
- Measurement cup
- Mixing cup
- Mixing spatula
- Channel wax (33 strips  $\varnothing$  7 mm and 30 strips  $\varnothing$  3 mm)
- Timer
- Flask filling basket

### Palajet supplementary unit

- Injection cylinder including cover
- 2 flasks (Palajet duoflask, Palajet singleflask)
- 12 injection cylinder cartridges with 12 cover cartridges and 50 sealing disks
- Timer
- Flask filling basket

Timer

## Palajet accessories

- Palajet singleflask
- Palajet duoflask
- Injection cylinder including cover
- Timer
- Flask filling basket
- Spacer

## 12. Name plate

<b>Heraeus Kulzer GmbH</b>		0600-07-2011
Grüner Weg 11, D-63450 Hanau		
P = 5-10 bar	P <sub>N</sub> = 4 bar	
66020450	Palajet®	
		CE

## 13. Contacts in the countries

LAND / COUNTRY	NAME / ADRESS
<b>Deutschland / Germany</b>	Heraeus Kulzer GmbH, Produktsupport Labor, Grüner Weg 11, 63450 Hanau Tel. / Phone +49 (0) 61 81 / 35-58 94 oder / or 35-47 73, Fax +49 (0) 61 81 / 35-59 93
<b>Australien / Australia</b>	Heraeus Kulzer Australia Pty. Ltd., Rydecorp, Unit 6, 2 Eden Park Drive, Macquarie Park NSW 2113 Tel. / Phone (02) 8422 6100, Fax (02) 9888 1460
<b>Brasilien / Brazil (America Sul / America del Sur / South America)</b>	Heraeus Kulzer South America Ltda., Rua Cenzo Sbrighi, 27 – sala 42 – São Paulo – SP – CEP 05036-010 Tel. / Phone +55 11 36650500, Fax +55 11 36650521
<b>China / China</b>	Heraeus Kulzer Dental Ltd., 1585 Gu Mei Road, 200233 Shanghai Tel. / Phone +86 21.649.58488, Fax +86 21.649.51732
<b>Frankreich / France</b>	Heraeus, Division Dentaire, 12, Avenue du Québec, Villebon - B.P.630, 91945 Courtaboeuf Cédex Tel. / Phone +33 169.18.48.85, Fax +33 169.28.78.22
<b>Großbritannien / United Kingdom</b>	Heraeus Kulzer Ltd., Albert Road, Northbrook Street, Newbury, Berkshire, RG14 1DL Tel. / Phone +44 1635.30-500, Fax +44 (0) 1635 524622
<b>Indien / India</b>	Heraeus Kulzer Dental India Private Ltd., Heraeus Kulzer House 344/2, Ladoo Sarai 110030 New Delhi Tel. / Phone +91 116512849, Fax +91 116512869
<b>Italien / Italy</b>	Heraeus Kulzer S.r.l., Via Console Flaminio 5/7, 20134 Milano Tel. / Phone +39 02210.09.41, Fax +39 02210.09.42-83
<b>Japan / Japan</b>	Heraeus Kulzer Japan Co., Ltd., TSK Bldg. 2F, 8-13 Hongo 4-chome Bunkyo-ku, Tokyo 113-0033 Tel. / Phone +81 35803.21-51, Fax +81 35803.21-50
<b>Mexiko / Mexico</b>	Heraeus Kulzer Mexico S.A. de C.V., Homero 527, 301 y 302 Co. Pol., 11560 Mexico Tel. / Phone +52 55.5531-5549, Fax +52 55.5255-1651
<b>Niederlande / Netherlands</b>	Heraeus Kulzer Benelux B.V., Postbus 986, NL-2003 RZ Haarlem Tel. / Phone +31 23.543.42-50, Fax +31 23.543.42-55
<b>Nordamerika / North America</b>	Heraeus Kulzer, LLC, Headquarters, 300 Heraeus Way, South Bend, IN 46614 Tel. / Phone 1-800-431-1745, Fax 1-800-522-1545
<b>Österreich / Austria</b>	Heraeus Kulzer Austria GmbH, Nordbahnstr. 36/2/4/ Top 4.5, A-1020 Wien Tel. / Phone +43 1.408.09.41, Fax +43 1.408.09.41-70
<b>Schweiz / Switzerland</b>	Heraeus Kulzer Schweiz AG, Ringstrasse 15A, CH-8600 Dübendorf Tel. / Phone +41 43.333.72-50, Fax +41 43.333.72-51
<b>Skandinavien / Scandinavia</b>	Heraeus Kulzer Nordic AB, Box 437, SE-191 24 Sollentuna, Hammarbacken 4B Tel. / Phone +46 8585.777.55, Fax +46 8623.14.13
<b>Spanien / Spain</b>	Heraeus S.A., Forjadores, 16, Prado del Espino, 28660 Boadilla del Monte, Madrid Tel. / Phone +34 91358.03-75, Fax +34 91358.03-68

## 14. Document history

07.2011 first edition